

Image-Pro[®]
by Media Cybernetics, Inc.

IP

Version 10

Image Analysis Starts with Image-Pro®

Image-Pro from Media Cybernetics® delivers the power and flexibility that you need in a modular image analysis platform, to address any application.

- Add Image-Pro Modules to expand the functionality
- Download APPS from the Image-Pro APP Store to simplify your analysis
- Customize the system to best meet your image analysis needs

Regardless of the application, Image-Pro functionality can be easily expanded or contracted to target your exact analysis need and deliver results.

To learn how Image-Pro can work for you, obtain a free 14-day trial of Image-Pro at: www.mediacy.com/freetrial



The **Image-Pro** platform offers simplified, flexible 2D solutions for any application.

Modules

Augment Image-Pro with additional modules to expand the product functionality

Apps

Simplify your analysis workflow by implementing a standard APP or develop your own solution.






Customization

Configure product interface toolbars, and reports to meet your specific needs.

Image-Pro Applications and Solutions

Go beyond the 2D analysis capabilities of Image-Pro with the addition of Image-Pro Modules. Modules open a new level of functionality for Image-Pro. Include them in your initial configuration or add them later, as your needs change.

To determine the best configuration, visit our application pages for suggested modules and other application specific information to optimize Image-Pro for your work.



Life Science


- Cell Counting
- Counting Markers within Cells
- Angiogenesis
- Wound Healing
- Cell Tracking
- 3D Cell Analysis
- 3D Colocalization

Manufacturing

- Inspection with Extended Depth of Field
- Quality Assurance
- Aerosol Spray Analysis
- Intermetallic Corrosion
- Particle Size Analysis
- 3D Visualization and Volume Analysis


Digital Pathology

- Percent Area Staining
- Cell Classification
- Microglia/Brain Tissue
- Color Deconvolution
- Whole Scanned
- Slides/Tissue Analysis




Materials Sciences

- Ceramics, Composites and Plastics
- Metals
- Cosmetics
- Glass, Films and Coatings
- R&D of New Materials
- Food Quality Assurance
- Fiber Analysis
- Particle and Porous Materials Analysis
- Layer Thickness




Natural Resources

- Entomology
- Forestry and Plant Disease
- Agriculture
- Remote Sensing
- Population Analysis
- Otolith Analysis and Morphology
- Egg and Larvae Counting
- Fisheries



Electron Microscopy

- Particle Size Analysis Metals
- Fiber Analysis
- Phase Analysis
- Image Tiling
- FIB-SEM 3D Volume Rendering and Analysis




Security

- Forensic Image Optimization
- Fiber Analysis
- Tool Marks Analysis
- Remote Sensing

Image-Pro Modules

2D Image Capture Module

The Image-Pro Capture Module goes far beyond image capture and camera support. It contains a variety of specialty capture options such as High Dynamic Range (HDR) capture, Live Tiling, Live Extended Depth of Field (EDF), Live Measurements, Live Image Comparison, and Live Calibrated Grids.



3D Module

Add the Image-Pro 3D Module to your platform and unleash a new level of image analysis for your research. Advanced 3D visualization and analysis options based on the same philosophy as the Image-Pro platform, allow you to perform complex analysis methods without learning a new platform.

License Options

Image-Pro offers two different licensing options to meet your organizational and laboratory needs. Available either as a single licensed copy or a network license. An Image-Pro single seat license is contained on the PC where Image-Pro has been installed, while in contrast, a network license supports multiple instances of Image-Pro across your internal network, through a license hosted on a central server.

	Single Seat License	Network License
License Type	USB Hard License (Dongle)	Soft License on Network Server
Base Functionality	Image-Pro	Image-Pro (per seat)
Number of Seats	1	>1
Modules Supported	Linked to USB License at Login	Assigned through a License Manager

Network licenses implement the modularity of Image-Pro among the instances running on the internal network. Each instance is an independent copy of the Image-Pro Software, each capable of having a unique configuration. Through the server licensing tool, the system can be reconfigured to move a module from one instance of Image-Pro to another or add a new module to a single instance or multiple instances.

Example of Image-Pro Network Configuration:




Image-Pro Modules Feature Chart

Image-Pro

Image Capture Module

3D Module

Object Tracking Tools

- Manual and automatic object tracking tools
- Make morphological measurements on the objects being tracked
- Correlation tracking
- Define reference track
- Compare tracks

-
-
-
-
-

3D Rendering and Analysis

- 3D filters
- Additional 3D filters (Top-hat Plus, Minus & Ring, Local Max & Min, & Gradient)
- 3D viewer
- Enhanced 3D viewer and volume rendering
- Clipping plane viewing
- Slicer viewing
- Iso-Surface based object analysis
- Volume "Pseudo-surface" based object analysis
- Guided Segmentation for complex samples
- Automated 3D volume measurements
- 3D object based manual measurements
- Batch processing of 3D volumes
- Simple Key Frame Animation movie making

-
-
-

-
-
-
-
-
-
-
-

Image Analysis Tools

- Data Collector collects data from multiple images from various tools
- Data Collector Class Statistics
- Surface Plot displays image intensity values in a 3D plot
- Measurement equation editor
- Export data to Microsoft® Excel
- Data graphing features include tools to dynamically link to image overlays
- Live data collecting and charting

-
-
-
-
-
-
-

Auditing and Authentication Features

- Audit trail
- Image and file signature "fingerprint" images for authenticity

-
-

Image Management and Reporting

- Annotation tools
- Create customized reports with images, data, and text
- Export measurements to statistical and spreadsheet package

-
-
-

Automation & Customization Tools

- Macro programming tools
- Graphical macro programming
- Customize workflow toolbars
- Customizable Apps
- Customize your application with the interactive menu editor
- Intelligent macro editor with command recognition
- Time-saving "Extras" macros in Macro Menu
- Downloadable Apps from the Apps Center
- Macro browser
- Built in Batch Processing Features
- Batch process - Folder Monitoring
- Link to online tutorials
- Customizable Ribbon

-
-
-
-
-
-
-
-
-
-
-
-
-

Operating System Support

- Windows® 7 Professional and up, (64-bit)
- Windows® 8/8.1 Professional and up, 64-bit
- Windows® 10 Anniversary Update, Professional and up, 64-bit
- Windows® 32-bit version are supported but not recommended

-
-
-
-

-
-
-
-

-
-
-

File Format Support

<http://www.mediacy.com/imagepro/system-requirements>

-



System Requirements

	Recommended	High Performance for 3D Module
Operating System	Windows 10 (64bit) - Anniversary Update, Professional and up Windows 8.1 (64bit) - Professional and up Windows 7 (64bit) - SP1, Professional and up *32bit versions are supported but not recommended.	Windows 10 (64bit) - Anniversary Update, Professional and up
Free Disk Space	Multi high-speed SATA hard disks or SSDs; 8GB free on installation drive + free space for images (20 GB+)	Multi high-speed SATA hard disks or SSDs configured in RAID 5; 16GB free on installation drive + free space for images (1 TB+)
Processor	2.8 GHz CPU Intel quad-core processor or better	3.0 GHz CPU Intel quad-core processor or better
RAM	16 GB memory or higher	48 GB memory or higher
Graphics Card	NVIDIA GeForce GTX Cards with 4GB graphics memory & Open GL 4.2 or higher	NVIDIA GeForce GTX Cards with 8GB graphics memory & Open GL 4.2 or higher
Drives	DVD-ROM drive if optional DVD purchased	DVD-ROM drive if optional DVD purchased
USB Port	Required for hard licenses and offline license activation	Required for hard licenses and offline license activation
Internet Connection	Internet connection required for online services*	Internet connection required for online services*
Internet Browser	Internet Explorer version 9 or higher	Internet Explorer version 9 or higher

Resources

Maintenance

Receive installation and configuration troubleshooting from our technical support engineers via phone or email. Purchase an Image-Pro Maintenance Agreement to receive a full year of software updates.

Training

Our image analysis training courses provide valuable information you can apply towards your imaging work. Whether you are new to imaging or are a more experienced imaging professional, our courses are designed to help you make optimum use of the latest imaging software techniques.

Get Started with Free Video Tutorials

As an Image-Pro user, you'll have access to a growing library of free 3 to 5 minutes video tutorials. Learn timesaving tips and familiarize yourself with the wide range of analysis tools available to you - without ever having to leave your computer.

Image-Pro forums

When you choose Image-Pro software, you instantly join a community of image analysis professionals. Join the free Image-Pro users forum to ask questions and learn from other users.

For more information, visit us at www.mediacy.com

